

Accuhaler™



Steps for using the Accuhaler™

1. Open the cover by pushing the thumb grip
2. Holding the device horizontally, load the dose by sliding the lever until it clicks
3. Breathe out gently, away from the device
4. Place lips around the mouthpiece without biting and form a good seal, keeping the inhaler horizontal.
5. Breathe in steadily and deeply
6. Hold breath for 5-10 seconds, or as long as comfortable
7. Breathe out gently, away from the device
8. Wipe the mouthpiece with a dry tissue
9. Close the lid by sliding the thumb grip until it clicks back into place.

If an extra dose is required, repeat steps 1 to 9.

When is the Accuhaler™ empty?

The Accuhaler™ is empty when the dose counter on the top shows a zero. The last five doses will appear in red.

How do I clean and store the Accuhaler™?

Wipe the inhaler with a clean, dry tissue after use. Do not wash any part of the Accuhaler™ or get it wet. Store in a dry place below 30°C, away from direct heat or sunlight.

Check the expiry date and ask your pharmacist about disposal of medicines.

Note: Refer to manufacturer's instructions for additional information on this inhaler device.