



Ellipta®



## Steps for using the Ellipta®

1. Slide the cover down until you hear a click
2. Breathe out gently, away from the device
3. Place lips around the mouthpiece without biting and form a good seal, keeping the vents clear
4. Breathe in steadily and deeply
5. Hold breath for 5-10 seconds, or as long as comfortable
6. Breathe out gently, away from the device
7. Wipe the mouthpiece with a dry tissue
8. Slide the cover to close the device

**Note: The device should not be shaken at any time**

## When is the Ellipta® empty?

The Ellipta® is empty when the counter shows a zero. Opening and closing the cover without inhaling results in loss of dose.

## How do I clean and store the Ellipta®?

Wipe the inhaler mouthpiece with a clean, dry tissue before closing the cover.

Keep in a cool dry place below 30°C.

**Check the expiry date and ask your pharmacist about disposal of medicines.**

Note: Refer to manufacturer's instructions for additional information on this inhaler device.