



## Steps for using the Ellipta

1. Check the dose counter to make sure the inhaler isn't empty.
2. Slide the cover down until you hear a click. The click means the dose is loaded and the inhaler is ready to use.
3. Breathe out gently, away from the inhaler.
4. Place lips around the mouthpiece without biting and form a good seal. Keep the vents clear.
5. Breathe in quickly and deeply.
6. Take the inhaler out of your mouth and hold breath for 5-10 seconds, or as long as comfortable.
7. Breathe out gently, away from the inhaler.
8. Wipe the mouthpiece with a dry tissue.
9. Slide the cover closed.

**Note: Do not shake the inhaler.**

## When is the Ellipta empty?

The inhaler is empty when the counter shows a zero. If you open and close the cover without inhaling the medicine, you lose the dose.

**Note: Dispose of the inhaler one month after you first open it or when the counter shows a zero, whichever comes first.**

## How do I clean and store the Ellipta?

Wipe the inhaler mouthpiece with a clean, dry tissue before closing the cover. Do not wash any part of the inhaler or get it wet. Keep in a cool dry place below 30°C, away from moisture, heat or sunlight. Do not store it in the bathroom, near a sink or a window, or in the car.

**Check the expiration date on your medicine.** Write the date the inhaler should be discarded on the label in the space provided. Ask your pharmacist how to safely get rid of any medicine you don't need.

**If you have any questions about using this inhaler, talk to your doctor, nurse or pharmacist.** Your inhaler will come with instructions. Always read these instructions for important information on how to use it.